English	Instructions for use: 3D Printed Surgical Guides	2 – 4
Deutsch	Gebrauchsanweisung: 3D gedruckte Bohrschablonen	5 – 7
Français	Mode d'emploi : Guides chirurgicaux obtenus par impression 3D	8 – 10
Italiano	Istruzioni per l'uso: Guide chirurgiche stampate in 3D	11 – 13
Español	Instrucciones de uso: Guías quirúrgicas impresas en 3D	14 – 16
Português	Instruções de utilização: Guias cirúrgicas impressas em 3D	17 – 19
Svenska	Bruksanvisning: 3D-skrivna kirurgiska guider	20 – 22
Nederlands	Gebruiksaanwijzing: 3D-geprinte boormallen	23 – 25
Čeština	Návod k použití: 3D tištěné chirurgické šablony	26 – 28
Русский	Инструкция по применению: хирургические направляющие, изготовленные методом 3D печати	29 – 32

Manufacturer / Hersteller / Fabricant / Produttore / Fabricante / Fabricante /
Fabricante / Tillverkare / Producent / Výrobce / Изготовитель

Createch Medical , Poligono Kurutz Gain P3B
20850 Mendaro / Spain , www.createchmedical.com

Caution: U.S. federal law restricts this device to sale by or on the order of a dental professional.

1. Product description

The Createch Medical 3D printed surgical guide is a patient-customized product made from bio-compatible acrylic. The 3D printed surgical guide is designed by the customer using a Straumann-approved 3D planning and design software (coDiagnostiX).

The surgical guide incorporates stainless steel, PEEK or titanium T-sleeves that define the position, orientation and depth of the drill used to create the surgical implant orifices.

2. Intended use

The 3D printed surgical guide is a 3D printed device that includes stainless steel, PEEK or titanium precision-made T-sleeve(s) for assisting inguided implant surgery. The intended use of the surgical guide is to aid clinicians in dental implant positioning after they have pre-planned the implant position and the position of the drill sleeves in a validated implant planning software.

3. Indications

Surgical dental implant guided surgery using the 3D printed surgical guide for locating the required drill position can be applied to defined cases identified by dental professionals, taking into consideration all the indications for implant treatment that apply in each patient situation. In all cases the anatomy, bone quality, guide support features and relationship to adjacent structures should always be taken into account during the pre-surgery planning.

Surgical dental implant guided surgery is indicated for assisting in implant placement in the following situations;

- single missing tooth,
- partially edentulous jaw or several missing teeth,
- edentulous jaws

4. Contraindications

The use of a surgical guide is contraindicated in the following circumstances:

- ☞ Patients not medically compliant with oral surgery,
- ☞ Patients whose oral situation does not fulfill the criteria for achieving a safe and functional outcome after surgery,
- ☞ Patients who are allergic or hypersensitive to stainless steel, PEEK, titanium or acrylic materials.

5. Side effects, interactions and precautions

Disclaimer

The 3D printed surgical guide has been tested and validated on average designs created using the approved 3D planning and design software coDiagnostiX (for details please refer to the coDiagnostiX Instructions for use). If the customer chooses to design a 3D printed surgical guide with dimensions outside the coDiagnostiX average range, Createch Medical cannot give any guarantee, and does not assume any liability for, the stability of the result.

Caution

This is a single-use product.

Warning

Use of a non-sterile device during surgery may lead to infection of tissues or infectious diseases.

6. Cleaning and sterilization

The surgical guides are made from polymerized acrylic materials.

When cleaning is required, we recommend cleaning the surgical guides only with approved materials (e.g. isopropyl alcohol (IPA) with less than 5 minutes submersion. To prevent the surgical guides from absorbing IPA residue, do not place them in a closed container or bag until the IPA evaporates completely).

The restoration should be placed in an accessory cassette and double-packed in common sterilization wraps (paper/film/bags). Steam sterilize according to the following parameters.

Caution: Ensure that any accessories used for wrapping or packaging are CE approved.

Material	Method	Conditions
Acrylic 3D printed surgical guide	Autoclave, moist heat Fractionated vacuum	North America 132 °C (270 °F), 4 min 30 minutes drying time
		Europe 132 °C (270 °F), 3 min 30 minutes drying time

Warnings

- ☞ Always keep the collar of the T-sleeve(s) uppermost in order to ensure that the sleeves remain stable during sterilization

- ☞ Support the guide with the wraps to prevent distortion during sterilization
- ☞ Do not place anything on top of the guide during sterilization
- ☞ After sterilization let the guide cool for 30 minutes before moving it

Please note:

If visible signs of moisture are present (damp spots on sterile packaging, pooled water in the load), please repackage and re-sterilize with a longer drying time.

Caution: Use devices immediately after sterilization. Do not store sterilized devices.

7. Procedure

Design of the 3D printed surgical guide – steps:

1. Before designing a 3D printed surgical guide, to the customer must have the patient's CT scan so that the bone and the tissue situation can be assessed and identified.
2. After assessing the patient's situation, the customer plans the position and the type of implant(s) with respect to the anatomy. The planning is done using coDiagnostiX software; the customer should be aware of all the limitations, warnings and advice stated in the coDiagnostiX IFU.
3. The 3D printed surgical guide is designed after the planning of the implant positioning. The design should take account of all the relevant patient-specific customized parameters. Windows need to be created to assess the planned placement of the 3D printed surgical guide and to use T-Sleeve(s) for the 3D printed surgical guide holes. The T-Sleeve(s) are inserted in the 3D printed surgical guide during production. The 3D printed surgical guide is delivered from Createch Medical to the receiving party with the following content:
 - 010.6121 T-sleeve guided height 6 mm, Ø 2.8 mm
 - 010.6122 T-sleeve guided height 6 mm, Ø 2.2 mm
 - 010.6123 T-sleeve guided height 5 mm, Ø 5 mm
 - 010.6124 Template Fixation Pin, Ø 1.3 mm
 - 010.6125 T-Sleeve for Template Fixation Pin
 - 010.6198 T-sleeves self-locking guided height 5 mm, Ø 5 mm
 - Drill T-Sleeves, T-Sleeve for Template Fixation Pins and Template Fixation Pins from Neodent® and MEDENTiKA® implant may also be delivered with the 3D printed surgical guide. For additional information about

the use of the Straumann® Group Neodent® and MEDENTiKA® dental implant systems, Drill T-Sleeves, T-Sleeve for Template Fixation Pins, Template Fixation Pins and implant products, please contact the relevant manufacturer's customer service department

- After the design is completed, the dental professional can send the design to the Createch Medical centralized production using the **caseXchange** option in the coDiagnostiX software.

Warnings

- Follow all the coDiagnostiX recommendations to produce an effective guide.
- Ensure that the coDiagnostiX drilling protocol is printed and sent to the dentist/surgeon before the surgery.
- Clinicians should be aware of potential areas for inaccuracies within the 3D printed guide and review the source images with the 3D printed guide to ensure that the guide has the planned geometry.
- Ensure that observation windows are made in the guide so that the dental surgeon can align the guide correctly before drilling.
- Ensure that T-sleeves are in place for each of the drilling positions.
- The guide designer should add a patient ID to the 3D printed surgical guide when designing with coDiagnostiX in order to prevent mix-ups.
- Avoid long distances between supporting anatomy and the drill sleeve.
- For drill T-sleeve situations that are supported on only one end, ensure that there is a sufficiently secure connection to the drill T-sleeve and that the supporting end can be securely held during drilling.

Dental surgeon - Recommendations prior to usage of the surgical guide:

- Keep the 3D printed surgical guide away from direct sunlight and/or elevated temperatures to prevent deformation.
- Verify the alignment of the T-sleeves with your implant plan. Do not use the 3D printed surgical guide if this check does not have a positive result.
- We recommend that all the materials used (3D printed surgical guide, T-sleeves) should be checked for damage. Use only materials in perfect condition.
- The 3D printed surgical guide plan and design cannot be modified by Createch Medical, but only printed as specified. Createch Medical will neither substitute nor replace a dental professional's plan and design for the surgical guide.
- Ensure that the coDiagnostiX drilling protocol is available and corresponds to the planned implants and patient.

Dental surgeon during surgery

Caution

- The dental professional should secure the 3D printed surgical guide and all its accessories against aspiration during trial of the guide intra-orally and during the drilling procedure.
- If the dental professional notices the detachment of any part from the 3D printed surgical guide, the procedure should be stopped and the part removed immediately in order to avoid potential swallowing.
- Sterilize the guide according to Section 6.
- Make sure that the surgical guide can be properly positioned and stabilized on the jaw without gaps during use

8. Further information

- Instructions for use: Non-sterile surgical instruments and prosthetic auxiliaries (701124/en)*
- Straumann Guided Surgery, Basic Information (152.753/en)*

9. MRI safety information

The Createch Medical 3D printed surgical guides and sleeves have not been evaluated for safety and compatibility in the magnetic resonance environment. The Createch Medical 3D printed surgical guides have not been tested for heating or migration in the magnetic resonance environment.

10. Please note

Practitioners must have knowledge of dental implantology and instruction in the handling of the Createch Medical product described herein ("Createch Product") for using the Createch Medical Product safely and properly in accordance with these instructions for use.

The Createch Product must be used in accordance with the instructions for use provided by the manufacturer. It is the practitioner's responsibility to use the device in accordance with these instructions for use and to determine whether the device matches the individual patient situation.

The Createch Product is part of an overall concept and must be used only in conjunction with the corresponding original components and instruments distributed by Institut Straumann AG, its ultimate parent company and all affiliates or subsidiaries of such parent company ("Straumann"). Use of products made by third parties, which are not distributed by Straumann, will void any warranty or other obligation, express or implied, of Straumann.

11. Validity

Upon publication of these instructions for use, all previous versions are superseded.

© Createch Medical, 2021. All rights reserved.







12. Availability


Some items of the Straumann® Dental Implant System are not available in all countries.

Symbols

The following table describes the Symbols that may be printed on the packaging label.

Please refer to the packaging label for the applicable symbols related to the product.

Symbol	Symbol Description	Symbol Source
	Consult instructions for use. Please follow the link to the e-IFU www.ifu.createchmedical.com	ISO 15223-1
	Manufacturer	ISO 15223-1
	Date of manufacture	ISO 15223-1
	CE marking is the manufacturer's declaration that the product meets the requirements of the applicable EC legislation	CE Certification Mark
	Catalogue number	ISO 15223-1
	Batch code	ISO 15223-1
	Serial number	ISO 15223-1
	Do not re-use	ISO 15223-1
	Non-sterile	ISO 15223-1
	Sterile	ISO 15223-1
	Sterilized using Ethylene Oxide	ISO 15223-1
	Sterilized using aseptic processing techniques	ISO 15223-1
	Sterilized using irradiation	ISO 15223-1
	Sterilized using steam or dry heat	ISO 15223-1
	Do not re-sterilize	ISO 15223-1
	Do not use if package is damaged	ISO 15223-1
	Keep away from sunlight	ISO 15223-1

Symbol	Symbol Description	Symbol Source
	Quantity symbol	Institut Straumann AG
	Use-by date	ISO 15223-1
	Temperature limitation (2°C-8°C/36°F-46°F)	ISO 15223-1
	Keep dry	ISO 15223-1
	Caution	ISO 15223-1

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