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DATA SHEET :**NOTE FOR THE LABORATORY
CrCo:Additive****Technical characteristics:**

Composition (average in mass - %):

Co – 63,9

Cr – 24,7

W – 5,4

Mo - 5

Density: 8,5 g/cm³

Expansion coefficient: (25°-500°) **14.1 µm/K**

Vickers hardness: 460HV10

Modulus of elasticity Gpa: 210

1. Polishing and brushing:

Brush the surface of the structure with Perlablast® micro, then polish it with a rubber polisher and brush it with a suitable polishing paste. Finally, clean the surface with steam.

2. Oxidation:

Oxidation is not recommended, if it has been evaluated that oxidation is necessary (perform 5 min at 900°C, in a vacuum environment).

3. Bonding:

Bonding is recommended in order to improve adhesion properties. The phase in which it should be applied shall depend on the brand of bonding to be used (see manufacturer's instructions).

4. Opaque:

Apply different opaque coatings (2-3 layers) with **approx. 10 min. of pre-drying.**

5. Ceramic coating:

Use ceramics with an adequate CET. Follow the instructions provided by the manufacturer of the ceramic. Slow cooling is recommended (cooling phase up to 600°C approx.). Avoid excessive thicknesses of the ceramic.

Acrilic coating:

The instructions of the coating manufacturers should be followed.

6. Polishing and brushing:

Polish the metal surface with a rubber polisher (without touching connections)

Polish the connections with polishing-pearl.

Guarantee:

The user is responsible for correct use and treatment. As a general rule, the instructions provided by the manufacturer of the ceramic should be followed. This document is for guidance purposes only.

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