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DATA SHEET:**NOTE FOR THE LABORATORY – 7
CoCr:PK****Technical characteristics:**

Composition (average in mass - %):

Co - 60

Cr - 29

W - 9

Si - 1,45

Fe - 0,15

Mn - 0,25

C - 0,01

Ni <0,10

Density: 8.3 g/cm³

Expansion coefficient: (20°-500°) **14.2 µm/K**

Elongation: 10%

Vickers hardness: 275HV10

Modulus of elasticity Gpa:245

Bonding:

According to the manufacturers of ceramics and in order to improve its adherence properties, bonding must be applied. The phase in which it should be applied shall depend on the brand of bonding to be used (see manufacturer's instructions).

PHASES:**1. Polishing and brushing:**

Polish and brush the CoCr structure with a dental burr in a single direction.

2. Preparation of the CoCr structure:

Sandblast the structure with aluminium oxide, grain size approx. 110 µm, at 2-3 bars of pressure.

3. Clean the CoCr anatomical bridge:

Clean the structure with steam.

4. Oxidation of the metal:

Oxidation baking (firing) is recommended: rapid heating to 960°C (increases of 100°/min), empty between 580°C and 949°C.

5. Sandblasting and green oxide removal:

Sandblast the surface with aluminium oxide, granulation approx. 110 µm and carefully clean the surface in order to eliminate the green oxide.

6. Clean the CoCr anatomical bridge:

Clean the structure with steam.

7. Opaque:

Apply different opaque coatings (2-3 layers) with **approx. 10 min. of pre-drying.**

8. Apply the dentine:

Once the baking (firing) temperature has been reached, a **slow cooling** of 10-15 minutes (depending on the size of the bridge) down to 600°C is advised. Following the first baking (firing) of the dentine, **do not steam** the bridge to clean it, but instead use soap and water only.

9. Glazed

Once the baking (firing) temperature has been reached, a **slow cooling** of 10-15 minutes (depending on the size of the bridge) down to 600°C is advised.

Guarantee:

The user is responsible for correct use and treatment. As a general rule, the instructions provided by the manufacturer of the ceramic should be followed. This document is for guidance purposes only.

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